


# *Pumani*

bCPAP


## 2016 Catalogue



HADLEIGH HEALTH  
TECHNOLOGIES

# Pumani bCPAP Package

The Pumani bCPAP Package comes with one Pumani bCPAP device and a 'starter kit' of accessories:



Products included in the Pumani bCPAP Package	Quantity Per Package	Part Number
1. Pumani bCPAP device	1	3-597
2. Bottle	1	-
3. Patient Tubing	3	-
4. Hat Clips	6	-
5. End Cap	3	-
6. Power Cord (UK plug)	1	-
7. Bottle Tubing	1	-
8. Variety Pack of Nasal Prongs (6 prongs, sizes 0-5)	1	-
9. Large Hat	1	-
10. Medium Hat	1	-
11. Small Hat	1	-
12. User Manual & Repair Manual (not shown)	1	-
13. Spare Parts Kit (not shown)	1	-


Additional Pumani Accessories  
Available Separately

---

# Patient Accessory Items


---

1. Nasal Prongs
2. Hats
3. Hat Clips
4. Patient Tubing
5. End Cap
6. Fish Elbows



# Spare Items


- 7. Bottle
- 8. Bottle Tubing\*
- 9. Filter\* (not shown)



\*1 of each of these items are already included in the Pumani Spare Parts Kit

# 1. Nasal Prongs

The Nasal Prongs come in six sizes\*, and can be ordered in packs of one size or a in variety pack




Product	Quantity Per Pack	Part Number
Pumani Nasal Prong Pack, Size 0	6	3-520
Pumani Nasal Prong Pack, Size 1	6	3-521
Pumani Nasal Prong Pack, Size 2	6	3-522
Pumani Nasal Prong Pack, Size 3	6	3-523
Pumani Nasal Prong Pack, Size 4	6	3-524
Pumani Nasal Prong Pack, Size 5	6	3-526
Pumani Nasal Prong Pack, Variety	6	3-527

\*See pg. 11 for a Nasal Prong Size Guide

# 2. Hats


## 3. Hat Clips

The Hats come in three sizes\*: Small, Medium, and Large



Product	Quantity Per Pack	Part Number
Hat Pack, Small	10	3-533
Hat Pack, Medium	10	3-534
Hat Pack, Large	10	3-535

The Hat Clips secure the Patient Tubing to the Hat



Product	Quantity Per Pack	Part Number
Hat Clip pack (one size)	10	3-536


\*See pg. 11 for a Hat Size Guide

# 4. Patient Tubing

# 5. End Caps

# 6. Fish Elbows


---



Product	Quantity Per Pack	Part Number
Patient Tubing Pack	1	3-539



Product	Quantity Per Pack	Part Number
End Cap Pack	10	3-532




Product	Quantity Per Pack	Part Number
Fish Elbow Pack	10	3-530




# 7. Bottle

# 8. Bottle Tubing

---



Product	Quantity Per Pack	Part Number
Bottle	1	3-538




Product	Quantity Per Pack	Part Number
Bottle Tubing	1	3-587

# 9. Filter

---

The filter is placed under the pump to keep dust from entering the pump, and may require replacement every 1-2 years

Product	Quantity Per Pack	Part Number
Filter	1	3-541



# Nasal Prong & Hat Size Guide

---

Patient Weight Range	Nasal Prong Size
Less than 1,000 grams	0
1,000 grams to 1,250 grams	1
1,250 grams to 2,000 grams	2
2,000 grams to 3,000 grams	3
3,000 grams to 4,000 grams	4
Over 4,000 grams	5

Patient Weight Range	Hat Size
Less than 1,500 grams	Small
1,500 grams to 3,000 grams	Medium
Over 3,000 grams	Large

# Suggested Accessory Order Quantities

Accessory Items	Part Number	Quantity per Pumani device per year*
Pumani Nasal Prong Pack, Variety**	3-527	2
Pumani Nasal Prong Pack, Size 0	3-520	2
Pumani Nasal Prong Pack, Size 1	3-521	2
Pumani Nasal Prong Pack, Size 2	3-522	2
Pumani Nasal Prong Pack, Size 3	3-523	2
Pumani Nasal Prong Pack, Size 4	3-524	2
Pumani Nasal Prong Pack, Size 5	3-526	2
Hat Pack, Small^	3-533	3
Hat Pack, Medium^	3-534	3
Hat Pack, Large^	3-535	3
Hat Clip Pack	3-536	3
Patient Tubing**	3-539	47
End Cap Pack	3-532	1
Fish Elbow Pack	3-530	1
Bottle	3-538	1
Bottle Tubing	3-587	1
Filter	3-541	1

\*These estimates are based on the assumption that 50 patients are treated with one Pumani bCPAP device per year. You may need to adjust these quantities if your facility's patient load is significantly higher or lower than this estimate.

\*\*These items are labeled as single-use only; HHT recommends that you follow your facility's standard cleaning procedures if you chose to reuse these items at your own risk.

^HHT recommends that the hats are single-use only, as the hats will stretch after use and will not stay on the patient's head as well.



For ordering information, please contact Hadleigh Health Technologies:

[info@hadleighhealthtechnologies.com](mailto:info@hadleighhealthtechnologies.com)

[www.hadleighhealthtechnologies.com](http://www.hadleighhealthtechnologies.com)

Hadleigh Health Technologies, LLC

30 Castro Ave

San Rafael, CA 94901, USA

+1 (415) 454-3005